

# Disposable Medical Middle Drape

## 一、 Product name, model and specification

Product name : Disposable Medical Middle Drape

Model : A type、 compound B type、 sewing B type、 compound C type、 sewing C type、 compound D type、 sewing D type

Specification : width 200~5000mm , length 200~5000mm

## 二、二、 Name, address, production address, contact information and after-sales service unit of the registered person/manufacturer

1、 Name: Hunan Zhenheyikang Medical Instrument Co.,LTD

2、 Address: No.6 Building Jingxiang Energy, No.55 Xianguang east road, Gaoxin district Xiangtan, Hunan

3、 Production address:NO.6 Building, NO.8 Building,5th floor Jingxiang Energy, No.55 Xianguang east road,Gaoxin district Xiangtan, Hunan

4、 Contact number : +86 17711627084

5、 After-sales service unit: Hunan Zhenheyikang Medical Instrument Co.,LTD

三、 Production Registration Certificate No: Xiangtan Medical Products Administration  
Production Registration NO 20140004

四、 Product Registration Certificate No、 Product Technical Requirements No:  
XTFDA Certified 20210061

## 五、 Product description, intended use

The product is made of nonwoven fabric and / or plastic film and or absorbent paper composite or sewing processing. The product is disposable, non-sterile.

For health care in hospital beds or examination beds.

## 六、 Use method : Open the bag can be used, directly to treat the wound

## 七、 Contraindications, precautions, warnings and reminders

Product contraindications: None

Damaged packaging prohibited use, expired products prohibited use

Disposable use, please destroy after use

Avoid high temperature, humidity, direct sunlight

## 八、八、 Interpretation of graphics, symbols, abbreviations and other contents used in medical device labels



一次性使用



参考使用说明书



包装破损切勿使用

## 九、 Installation and use instructions or illustrations

This product does not need to be installed, open the bag can be used.

## 十、 Product maintenance and maintenance methods, special storage, transportation conditions, methods

Products should be stored in a non-corrosive gas and well-ventilated room, transported by ordinary transportation means, keep clean and pollution-free, rain and snow immersion.

## 十一、 Date of production, shelf life or expiration date

Valid for three years, production batch number, production date, expiry date please refer to the small package or qualification label.

## 十二、 Accessories list, including accessories, appendix, lossitage replacement cycle and replacement method description

This product has no accessories

## 十三、 Date and version of preparation or revision of the specification

This specification was revised in May 2023, version 2023.